

DGX 198

Sprague Rappaport stethoscope



DGX198 is a Sprague Rappaport stethoscope with two heads particularly versatile. The included equipment, allows to adapt the device to neonatal, pediatric and adult patients and makes it suitable for particular types of auscultations depending on the used configuration.

Specific features

- Extremely light
- Two-tubes structure
- Metal clips limit the friction of the tubes
- Equipped with adult, pediatric and infant bells
- Spare ear-tips (2 big and 2 small)
- Large and small spare diaphragms

Technical data

Large head diameter	45 mm
Small head diameter	31 mm
Tubes length	56 cm
Weight	235 ± 10g

Class I Med.Device according to Dir. 93/42/CEE

DGX 198 BLACK SPRAGUE
RAPPAPORT TYPE BLACK
CHROMED
DG21198C



CND Classification	C9005
Registration number	90164

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it