

## DGX 110

### Black basic phonendoscope



Phonendoscope simple, light and ultra-thin, it has a particular diaphragm suitable for auscultation during blood pressure measurement.

#### Specific features

Steel head  
Adjustable binaural arc

#### Technical data

Head diameter	47 mm
Weight	80 g
Material	Steel , PVC

Class I Med.Device according to Dir. 93/42/CEE

DGX 110 - NURSE SCOPE -  
PHONENDOSCOPE

**DG02110B**

**CND Classification** C9005

**Registration number** 90167



Rev.0 (25/07/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)