

DGX 110

Black basic phonendoscope



Phonendoscope simple, light and ultra-thin, it has a particular diaphragm suitable for auscultation during blood pressure measurement.

Specific features

Steel head
Adjustable binaural arc

Technical data

Head diameter	47 mm
Weight	80 g
Material	Steel , PVC

Class I Med.Device according to Dir. 93/42/CEE

DGX 110 - NURSE SCOPE -
PHONENDOSCOPE

DG02110B



CND Classification C9005

Registration number 90167

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it