

DGX 280

Pediatric Stethoscope



Stethoscope intended for chest auscultation of pediatric patients.

Specific features

Extremely light
Bright color to improve the approach to the child

Technical data

Large head diameter	37 mm
Small head diameter	26 mm
Weight	95 g

Class I Med.Device according to Dir. 93/42/CEE

DGX 280 PEDIATRIC
STETHOSCOPE

DG02280D



CND Classification C9005

Registration number 90171

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it