

# DGX 210

## Electronic thermometer



Electronic device for the determination of body temperature.

### Specific features

- Easy to use
- Acoustic signal for the end of reading
- Memory and fever warning
- Auto power off after 8 minutes
- Rigid protective case

### Technical data

Measuring range	32°C – 43°C
Precision from 35°C to 39°C	± 0,10°C
Display resolution	< 0,1 °C
Batteries	AG3 Alkaline or LR41 silver oxide
Batteries life span	2 years or 2500 uses
Dimensions	125x19x11 mm
Weight (batteries included)	9,8 g

Class IIa Med. Dev. according to Dir. 93/42/CEE

DGX 210 - ELECTRONIC  
THERMOMETER  
DG10210E



CND Classification

V03010102

Distributed by Spencer Italia s.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)