

DGX 091

Pen Light



Light for the inspection of eye and oropharyngeal cavity.

Specific characteristics

Stainless steel body
 Incandescent light
 Requires 2 batteries AAA not included

Technical data

Length	139 mm
Weight without batteries	13 g

Class I Med.Device according to Dir. 93/42/CEE

DGX 091 – STAINLESS STEEL
 PEN LIGHT WITHOUT
 BATTERIES



DG08091D

CND Classification	V9099
Registration number	797007

Rev Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it