

## LED LIGHT PEN



Light pen to be used by professionals. Thanks to the wide beam, constitutes a valuable aid for rescuers operating in conditions of poor visibility.

#### Specific Features

- Requires 3 batteries LR41 included
- Led light on the upper side of the pen
- Rigid case included
- Power button on the upper end to facilitate lighting
- Ergonomic and non slipping grip

#### Technical data

Dimensions	ø15x130 mm
Pen weight	26 g
Kit weight	41 g

LED LIGHT PEN - LUMINOUS PEN WITH  
LED  
DG08080D

Rev.0 (25/09/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italy Tel. +39.0521.541111 Fax +39.0521.541222  
[www.spencer.it](http://www.spencer.it)