

POP

Oxygen module for wall mounting



Oxygen module for wall recessed mounting, to be installed inside a medical vehicle on dedicated surfaces.

Specific features

Simple and safe use guaranteed by the adoption of gas-specific connection and quick release mechanism
Smooth surfaces for an easy cleaning
The locking mechanism, allows easy installation

Technical data

Overall dimensions	100x100x115 ±5 mm
Dimensions of the non-recessed part	100x100x35 ±5 mm
Connection standard	UNI
Materials	Nylon, Al, Brass

Compliant with

UNI 9507
UNI EN ISO 9170-1
EN 1789

Class IIb Med.Device according to Dir. 93/42/CEE

POP MODULE - OXYGEN
CONNECTION UNI 9507
OX04300A



CND Classification Z120309
Registration number 335394

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it