

POP

Vacuum module for wall mounting



Vacuum module for wall recessed mounting, to be installed inside a medical vehicle on dedicated surfaces.

Specific Characteristics

- Simple and safe use guaranteed by the adoption of gas-specific connection and quick release mechanism
- Smooth surfaces for an easy cleaning
- The locking mechanism, allows easy installation

Technical data

Overall dimensions	100x115x100 ±5 mm
Dimensions of the non-recessed part	100x35x100 ±5 mm
Connection standard	UNI
Materials	Nylon, Al, Brass

Compliant with

UNI 9507
UNI EN ISO 9170-1
EN 1789

Class IIb Med.Device according to Dir. 93/42/CEE

POP MODULE - VACUUM
CONNECTION UNI 9507

OX04301A



CND Classification Z120309

Registration number 416458

Rev.0 (26/03/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it