

POP

Three oxygen unit module for wall mounting



System for wall recessed mounting, with three medical oxygen output units, to be installed inside a medical vehicle on dedicated surfaces.

Specific Characteristics

- Simple and safe use guaranteed by the adoption of gas-specific connections and quick release mechanism
- Smooth surfaces for an easy cleaning
- The locking mechanism, allows easy installation

Technical data

Overall dimensions	300x115x100 ±5 mm
Dimensions of the non-recessed part	300x35x100 ±5 mm
Connection standard	UNI
Materials	Nylon, Al, Brass

Compliant with

UNI 9507
UNI EN ISO 9170-1
EN 1789

Class IIb Med.Device according to Dir. 93/42/CEE

POP MODULE 3 OXYGEN
SOCKETS UNI 9507
OX04309A



CND Classification Z120309
Registration number 335396

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it