

OXX 50

Vacuum probe



Vacuum probe for the connection of suction devices to a gas- specific socket compliant with the standard UNI, installed in systems for the distribution of medical gases.

Specific Characteristics

Allows quick connection and disconnection of devices from the sockets
Color coded for easy identification

Technical data

Dimensions	ø25x60 mm
Internal diameter	6 mm
Connection	F ¼"
Material	Nickel-plated brass

Compliant with

UNI 9507
EN ISO 9170-1

Class IIb Med.Device according to Dir. 93/42/CEE

VACUUM CONNECTION PROBE
UNI9507

OX05004A

CND Classification Z120309

Registration number 310705



Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it