

## OXX 50

Oxygen inlet



They represent the mobile male component, not interchangeable, designed to be inserted and locked into the outlet, ready to permit the suction of gas.

### TECHNICAL DATA

Material	Steel
Size	Ø24x 60mm
Weight	96 g
Gas	Medical oxygen
Flow rate	According to UNI EN ISO 7396-1
Terminal Unit	According to UNI 9507 and UNI EN ISO 9170
Latex free	

OX05013A                      OXX 50  
Oxygen inlet

CE 0123

UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>