

OXX 60

Oxygen inlet



They represent the mobile male component, not interchangeable, designed to be inserted and locked into the outlet, ready to permit the suction of gas.

TECHNICAL DATA

Material	Steel
Size	Ø24x 86mm
Weight	113 g
Gas	Medical oxygen
Flow rate	According to UNI EN ISO 7396-1
Terminal Unit	According to UNI 9507 and UNI EN ISO 9170 with hose holder ø7 mm
Latex free	
OX06000A	OXX 60 Oxygen inlet

CE 0123

UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>