

## O-Tube

### Oxygen connection tube



Tube intended for pneumatic connection of devices requiring the use of pressurized oxygen.

#### Specific features

Colored and marked for use with oxygen  
 Ideal for low pressure applications (max 14bar)  
 DEHP FREE

#### Technical data

Internal diameter	6,7 mm
External diameter	12,7 mm
Length	100 cm
Material	fabric reinforced PVC
Burst pressure	56 bar/23°C
Output connection	Quick coupling socket
Input connection	Threaded M1/4" G

DRIVING OXYGEN HOSE CM.  
 100  
**EV30020A**

Rev.0 (15/03/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. a socio unico Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)