

OXX 200

Compressor / tank



A compressor/tank group with the 12 V power source ready for use on the emergency vehicles.

In a small space the electro compressor, a non rebreathing valve, an electro valve to compensate the pressure and the pressure switch are wired together to guarantee the constant flow of compressed air.

Easy to install, it is supplied with power connections and connecting tube for the compressed air.

To obtain medical air it is essential that the correct filters and dryers are also assembled.

TECHNICAL DATA

Dimensions (mm)	315x190xh340
Weight (kg)	9,5
Power	12 Vcc – 16A
Absorption	192 W
Flow	70 L/min Rpm 2800
Use	30min
Max. pressure	3,5 bar
Tank Volume	4 Lt
Isolation Class	B
Protection grade	IP20
Cable Length	2 mt
Power switch	380 Vmax/60Hz 15 A max , 6bar
Limited use period	S2 30/5

OX02200A

OXX 200

