

## REK PRO

Clamp with thread for rail bar



Clamp to be used for fixation on rail bars of devices like sphygmomanometers, IV pole, flowmeters and other devices used inside the medical vehicle.

### Specific features

- Extremely compact and robust
- Quick coupling systems both for accessories and for the clamp itself
- High resistance to involuntary sideways movements thanks to the use of inserts with high adhesion at the points of contact with the rail bar
- Can be used with rail bars complying with the standard UNI EN 19054

### Technical data

Dimensions	67x60x90 mm
Materials	Nylon, Al
Weight	200 ± 5 g

Class I Med.Device according to Dir. 93/42/CEE

REK PRO THREADED CLAMP

OX07101B



CND Classification V9099

Registration number 446772

Rev.0 (25/07/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)