

REK PRO

Threaded clip for group bar G-Rail



Our highly qualified research and development team today is ready to present the fixation system for medical applications of tomorrow. Thanks to the use of advanced graphic and modelling systems like CAD software, our engineers have developed efficient and innovative solutions for any problem related to quick fixation. In fact, Rek is versatile, safe and durable and apart from the formal quality, it presents a second level of coherence: the technological and technical content.

Rek Pro allows to fix, organise and attach up various users, offering high stability even with vehicle in motion.

Using a ¼ " gas thread and the rubber connector diameter 6,5 (medical gas) or 8,5 (aspiration) it is possible to connect:

- Flowmeters
- Humidifiers
- Venturi aspirators
- Vacuum containers

In conformity with the norms:

UNI EN 1789 Emergency vehicles, their equipment and medical devices

UNI EN 19054 - Rail systems for supporting medical equipment

TECHNICAL DATA

Dimensions	67x60xh90
Weight	0,2 kg
Materials	Nylon - Aluminium

OX07102B

REK PRO 10G

