

MONITOX

Luminous Indicator for oxygen cylinders



When it is necessary to see the contents of the cylinders positioned in an area to which the operator does not have access (for example the cylinders are positioned in the driver's cab), Monitox is the ideal solution.

At any time you can see exactly how many bars remain in each of the cylinders.

The contents of the cylinder can be seen thanks to a green LED scale with red LED to indicate reserve.

The panel has a switch to block the flow and reduce use of the battery when the ambulance is not in use.

The reduced dimensions of the panel and the installation with Velcro® make it easy to install on any surface. The electronic manometer will allow you to view simultaneously the analogical pressure value of both scales. The dedicated cable connections are coded to avoid any mistake in assembly. Monitox is compatible with the automatic cylinder change Self Switch.

The electronic manometer (1/4 " gas) can be installed on good quality pressure reducers and on balanced pressure reducers PRX 05 Spencer.

Monitox is composed of an easy assembly kit which includes:

- Luminous panel with graduated printing
- Two electronic manometers
- Cables to connect panels and manometers
- 1 mt. Power cable

They are available also to assemble on personalised panels that can be connected to electronic manometers. This last solution is used also when one or more cylinders has to be monitored.

TECHNICAL DATA

Dimensions (mm)	82x70xh24
Weight (kg)	0,2
Power cable Length	1 mt
Length of sensor cable	2 mt
Power	12 Vcc
Absorbed power	0,5 A
Oxygen Manometer	0-300 bar – 1,6 C.P.
Manometer attachment	G 1/4"
Manometer diameter	50mm

OX07110A

MONITOX

