

OX-T

Oxygen tanks



The compressed oxygen tanks are the lightest of all steel type tanks. They have a pressure of 200 atm and are available in a wide range of sizes.

Medical gas oxygen tanks are supplied with

- chromed valve for oxygen
- steel disk for medical use
- Tank conformity certificate

TECHNICAL DATA

Material	Steel
Tank diameter	110 mm
Protection diameter	140 mm
Height	545 mm
Weight	4,36 kg
Liters	3 L
Oxygen volume at 200bar	600 L
Gas	Medical oxygen
Connection	UNI
Latex free	

OX10030E

OXYGEN TANK O₂ 3lt UNI
WITH FERRULE CAGE CUP
AND VALVE



The information contained in this technical document is compliant with the manufacturers declaration.

UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>