

OX-T

Oxygen tanks



The compressed oxygen tanks are the lightest of all steel type tanks. They have a pressure of 200 atm and are available in a wide range of sizes.

Medical gas oxygen tanks are supplied with

- chromed valve for oxygen
- steel disk for medical use
- cumulative test certificate
- charging connection UNI

TECHNICAL DATA

Material	Steel
Diameter	170mm
Height	920 mm \pm 5%
Liters	14 lt
Oxygen volume at 200bar	2800 lt
Gas	Oxygen
Latex free	

OX10140E

OXYGEN BOTTLE 14LT.
UNI WITH VALVE
AND PROTECTION



The information contained in this technical document is compliant with the manufacturers declaration.

UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>