

## OX-VA

Oxygen tank with integrated pressure reducer and flowmeter



Oxygen tank with flowmeter to be used in oxygen therapy.

### Specific features

- 1 liter oxygen tank with integrated flowmeter
- Can contain pressurized O<sub>2</sub> up to 200 bar
- High autonomy with reduced dimensions
- Flow adjustment from 0 to 15 l/min

### Technical data

Dimensions	ø140x360 ± 5 mm
O <sub>2</sub> Volume at 200bar	200 l
Connection standard	UNI
Output pressure	0 ÷ 2,8 bar
Output	Hose holder ø6 mm
Flow adjustment	From 0 to 15 l/min
Flow precision	±10% or ± 0,5 l/min
Manometer	0÷315 bar, cl 2,5, ø40 mm
Materials	Steel, PE, Chromed brass, Al, PC
Weight	2,5 ± 0,3 kg

Device marked CE according to dir. 93/42/CE

OXYGEN TANK 1 LT. W/VALVE

REDUCER AND PROTECTION

OX10012E

Distributed by Spencer Italia s.r.l.

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)