

OX-VA

Oxygen tanks with integrated
reducer valve and pressure
reducer/regulator



The compressed oxygen tanks are the lightest of all steel type tanks. They have a pressure of 200 atm, are available in a wide range of sizes and are complete of valve system, pressure reducer and dosage system.

Medical gas oxygen tanks are supplied with

- chromed valve for oxygen
- steel disk for medical use
- Tank conformity certificate

The information contained in this technical document is compliant with the manufacturers declaration.

TECHNICAL DATA

Material	Steel
Tank diameter	110 mm
Protector diameter	140 mm
Height	575 mm
Weight	5,16 kg
Liters	3 L
Oxygen volume at 200bar	600 L
Adjustable Flow	0-15 L/min
Outlet pressure	0÷3.5 bar
Gas	Medical oxygen
Connection	UNI
Latex free	

OX10031E

OXYGEN TANK 3lt UNI
W/VALVE.AND PRESSURE
REDUCER WITH PROT.

0123

UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>