

OX-Gr

Oxygen tank with pressure reducer

Flowmeter and humidifier



Oxygen tank complete of pressure reducer, flowmeter and humidifier intended for the treatment of respiratory failure.

Specific Characteristics

3 liters oxygen tank with integrated protection and large base to increase stability
 Can contain pressurized O₂ up to 200 bar
 High autonomy in reduced dimensions
 Piston pressure reducer with sintered filter
 Humidifier with high transparency to ensure an easy reading of the level

Technical data

Tank dimensions	ø110x545 ± 5 mm
Protection diameter	140 ± 5 mm
Assembly dimensions	595x145x260 ± 10mm
O ₂ volume at 200bar	600 l
Reducer inlet	UNI
Reducer output pressure	3,5/4 bar
Flowmeter regulation	From 0 to 15 l/min
Humidifier indicator of maximum	130 ml
Humidifier output	Hose ø7mm
Materials	Steel, PE, Chromed brass, Al, PC
Overall weight	5,2 ± 0,3 kg

Device subject to CE marking according to the dir. 93/42/CE

OX-GR OXYGEN BOTTLE 3L
 REDUCER/FLOWMETER/HUMIDIFIER

OX10032D

Rev.3 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it