

OX-Gr

Oxygen tank with pressure reducer
Flowmeter and humidifier



Oxygen tank complete of pressure reducer, flowmeter and humidifier intended for the treatment of respiratory failure.

Specific Characteristics

- 14 liters oxygen tank with integrated protection and large base to increase stability
- Can contain pressurized O₂ up to 200 bar
- High autonomy in reduced dimensions
- Piston pressure reducer with sintered filter
- Humidifier with high transparency to ensure an easy reading of the level

Technical data

Tank dimensions	ø170x920 ± 5 mm
Protection diameter	140 ± 5 mm
Assembly dimensions	920x170x260 ± 10mm
O ₂ volume at 200bar	2500 l
Reducer inlet	UNI
Reducer output pressure	3,5/4 bar
Flowmeter regulation	From 0 to 15 l/min
Humidifier indicator of maximum	130 ml
Humidifier output	Hose ø7mm
Materials	Steel, PE, Chromed brass, Al, PC
Overall weight	17,8 ± 0,3 kg

Device subject to CE marking according to the dir. 93/42/CE

OX-GR OXYGEN BOTTLE 14L
REDUCER/FLOWMETER/HUMIDIFIER
OX10142D

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it