

Oxygen tank with flowmeter and Oxygen socket



Indicative image



Oxygen tank equipped with flowmeter and oxygen socket.

Specific features

7 liters oxygen tank with integrated protection and base to increase stability
Can contain O₂ up to 200 bar

Technical data

O ₂ Volume at 200bar	1400 l
Input	UNI (W 21,7x1/14") Type B
Socket (output)	UNI
Socket output flow	140 ± 10 l/min
Output pressure	3,5 ± 0,5 bar
Flowmeter adjustments	0-1-1.5-2-3-4-6-9-12-15 l/min
Flowmeter output	Hose holder ø7mm
Manometer scale	0-300 bar
Materials	Steel, PE, Chromed brass, Al, PC

Device subject to CE marking according to dir. 93/42/CE

B-OX WITH 7L OXYGEN BOTTLE
AND PRESSURE REDUCER

OX10073B

Rev.0 (23/01/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it