

## PRX 05

### Membrane pressure reducer



The image can differ from the actual product



Pressure reducer with fixed output pressure suitable for direct connection to the oxygen tank. The pressure gauge allow to read the remaining pressure inside the oxygen tank.

#### Specific features

Input connection standard DIN  
 Safety valve  
 Output hose holder  $\varnothing 7\text{mm}$   
 Sintered bronze filter integrated in the device  
 Characterized by high and stable flow rate

#### Technical data

Maximum input pressure	200 bar
Output pressure	Factory set max 4 bar
Materials	Chromed brass, Steel, Bronze, Vytan
Manometer scale	0 – 300 bar precision class 2,5

#### Reference standard

ISO 10524-1:2006

Class IIb Med.Device according to Dir. 93/42/CEE  
 PRESS.RED.B3,5-150l/m-DIN  
 HOSE HOLDER.  $\varnothing 7\text{mm}$

**PR00350B**

**CND Classification**



Z120309

Distributed by Spencer Italia S.r.l

Rev.0 (05/04/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)