

PRX 08

Pressure regulator with dispenser



A professional pressure reducer with dosage selector. Manufacturing precision guarantees a perfect distribution of the flow supplied even when the tank is in a horizontal position. The necessity to operate with precision and safety will be answered even when operating conditions are not optimal. The protection of the colour-graded gauge, the regulation grip handles, the design intended to reach the highest ergonomic levels are translated into dynamism and maximum exploit of use.

TECHNICAL DATA

- Brass nickel structure
- Viton membrane and washer
- Sintered bronze filter 60 micron
- Graded gauge 40 mm diameter
- Proportioning device in non-slip-grip ABS
- Incoming pressure: 200 bar
- Outgoing pressure: 3,5 - 4 bar
- Outgoing precision: +/- 15% (from 20 to 200 bar)
- Pre-fixed safety valve to 7 bar
- Adjustable capacity with intervals of 20 bar
- UNI connection
- Outgoing connection complete with 6 mm diameter inlet
- Weight: 750 g

PR00308B PRESSURE REDUCER W/FLOW
SELECT 15l/m W/COVER
GAUGE - SHORT

CE 0434

UNCHECKED COPY - further revisions will be available on link <http://support.spencer.it/en/pdf>

Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it

Rev.0 (09/08/2013)

To verify the correctness of this document, please contact our Commercial Dept. at info@spencer.it