

PRX 08

Pressure reducer with flowmeter



Pressure reducer for direct connection to the oxygen tank. The pressure gauge allows easy reading of the pressure inside the oxygen tank.

Specific features

Input connection standard NF
 Safety valve
 Output with hose holder $\varnothing 6\text{mm}$
 Sintered bronze filter integrated in the device
 Flowmeter with steps with ABS knob

Technical data

Dimensions	100x110x130 max ± 5 mm
Weight	0,78 \pm 0,2 kg
Max input pressure	200 bar
Output pressure	Precalibrated max 4 bar
Flow adjustment (l/min)	0-1-1,5-2-3-4-6-9-12-15
Materials	Chromed brass, Steel, Bronze, Vyton
Manometer	0 – 300 bar precision class 2,5

Compliant with

ISO 10524-1:2006

Class IIb Med.Device according to Dir. 93/42/CEE

PRES.REG.DISP.15L/MIN NF
 W/GAUGE GUARD OUTLET
 $\varnothing 6\text{MM}$



PR00316B

CND Classification Z120309

Registration number 84055

Distributed by Spencer Italia S.r.l

Rev.0 (28/10/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it