

## FL-X

### Flowmeter/humidifier group



System including humidifier and a flowmeter intended to be connected to oxygen distribution systems.

#### Specific features

- Visualization of the gas flow delivered to the patient
- The humidifier helps to reduce the dryness of mucosas during administration of oxygen
- Humidifier with hose output for the connection of devices used for oxygen administration

#### Compliant with

- UNI EN 15002
- UNI EN ISO 8185

#### Technical data

##### Flowmeter with probe

Dimensions	85x150x35 mm
Nominal pressure	350 kPa
Flow	From 0 to 15 l/min
Deliverable gas	Oxygen
Materials	PC, anodized Al, Brass

##### Humidifier

Dimensions	ø 70x135 mm
Max. volume	130 ml
Output hose	ø 7 mm
Operating pressure	350 kPa
Max pressure	500 kPa
Max flow	20 L/min
Materials	PC, Al
Overall weight	250 g

Medical device CE marked according to dir. 93/42/CE

FL-X FLOWMETER /  
HUMIDIFIER GROUP

**FL08301B**

**Registration number** 21663

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)