

OXYMASK 50

Kit 10 medium concentration oxygen masks



Disposable face masks for oxygen therapy.

Specific Characteristics

Adjustable elastic
Semi transparent shell to improve patient monitoring
Anti-crush tube
Flexible aluminum nose clip to increase stability

Latex Free

Technical data

| | |
|----------------------------------|--------------------|
| Dimensions | 135x100x60 ± 10 mm |
| Material | PVC |
| Weight ^(1mask + tube) | 70 g |
| Tube length | 2 m |

Class IIa Med.Device according to Dir. 93/42/CEE

OXYMASK 50 MEDIUM
CONCENTRATION ADULT, 10 PCS
OX20139A



| | |
|----------------------------|-----------|
| CND Classification | R03010201 |
| Registration number | 104805 |

Rev.1 (09/09/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it