

SAM

Bag with oxygen bottle



Bag equipped with oxygen tank and other devices commonly used in ventilation procedures.

Included Equipment		Technical data	
RM30300A	COLLIN TONGUE FORCEPS	Dimensions	650x250x240 ± 20 mm
RM30700A	MOUNDER JAW SPREADER	Material	PVC, Nylon, Steel
RM10891A	CUFFED MASK SIZE 2	Weight empty	2,3 ± 0,1 kg
RM10892A	CUFFED MASK SIZE 3	Oxygen tank volume	2l
RM10894A	CUFFED MASK SIZE 4		
OX10022E	2L OXYGEN TANK W/FLOWMETER		
OX20120A	OXYMASK 50 MED. CONC. AD.		
OX20103A	OXYMASK 50 MED. CONC. PED		
OX20130A	ANTI-CRUSH OXYGEN TUBE		
RM10325A	GUEDEL DISPOSABLE CANNULA 40mm		
RM10320A	GUEDEL DISPOSABLE CANNULA 50mm		
RM10321A	GUEDEL DISPOSABLE CANNULA 60mm		
RM10322A	GUEDEL DISPOSABLE CANNULA 70mm		
RM10323A	GUEDEL DISPOSABLE CANNULA 80mm		
RM10324A	GUEDEL DISPOSABLE CANNULA 90mm		
BG02007A	ADULT PVC RESUSCITATOR BAG		
		Medical devices CE marked according to dir. 93/42/CE	
		SAM COMPLETE - BAG FOR OXYGEN BOTTLE	
		CB09014A	

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it