

## OXYPACK 1

Oxygen tank with pressure reducer

And manual resuscitator



Oxygen tank complete of pressure reducer and manual resuscitator.

### Specific Features

5 liters oxygen tank
Can contain pressurized O <sub>2</sub> up to 200 bar
High autonomy in reduced dimensions
Easy to carry thanks to the shoulder strap
Wide pocket with zip closure to store accessories
Hole with metal ring inside the pocket allows to pass an oxygen tube and to connect it to an oxygen mask
Pressure reducer with pressure gauge 0-300 bar with sintered filter and integrated flowmeter
Supplied with adult rubber disposable bag and oxygen tube for connection to the tank
PVC disposable mask size 3

### Technical data

Bag Dimensions <sup>(without tank)</sup>	640x270 ± 20mm
Side pocket dimensions	370x125x80 ± 5 mm
Bag weight	600 ± 10 g
Tank dimensions	Ø140x565 mm
Tank weight	8,3 ± 0,1 kg
O <sub>2</sub> volume at 200bar	1000 l
Connection standard	UNI
Output pressure	3,5/4 bar
Adjustable Oxygen Flow	0- 1-1,5-2-3-4-6-9-12-15 l/min
Pressure reducer weight	720 ± 10g
Pressure reducer output	Hose ø7mm
Resuscitation bag volume	1600 ml

### Med.Dev. CE marked according to the dir. 93/42/CE

OXYPACK 1 WITH OXYGEN  
CYLINDER 5 L PRESSURE  
REDUCER MANUAL RES. AND  
MASK  
**OX10206B**

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)