

AMBUJET P 800 D 20G

Portable suction device
with 20G wall fastener with auto recharge



Space-saving suction device to be used to remove fluids from the airway or respiratory support system. The device creates a negative pressure that draws fluids through disposable tubing that is connected to the collection container. The base performs as a fastener in accordance to EN 1789 standard.

Specific Characteristics

Low weight and size
Easy view of the vacuum set thanks to the led indicators
Classification ISO High Vacuum/ High Flow
Suction device protected from bumps and falls thanks to the protective shell
Mechanism of engagement and release from the base extremely rapid
Direct connection between canister and vacuum pump
Recharging of the battery is done simply by placing the device on the mounting system

Standard Equipment

10g fastener with quick unlocking system
Battery
Disposable canister
Disposable patient tube 1,8 m
Shoulder strap for easy transport

Optional accessories/Consumables

| | |
|----------|---|
| SC75012E | Lateral bags for Jet |
| SC75014D | Supply cable for cig.lighter 12vcc |
| SC73017E | 220V Power supply |
| RISC005E | Antibacterial filter |
| SC73016E | Sterile suction catheter 1,8m |
| SC73039E | 800 ml disposable canisters (10pcs set) |

Compliant with

ISO 10079-1
EN 1789

Technical data

| | |
|------------------------------------|---|
| Jet Dimensions | 200x225x285 mm |
| Jet Dimensions w/fastener | 200x250x320 mm |
| Suction device weight | 2,4 ± 0,1 kg |
| Fastener weight | 0,8 ± 0,1 kg |
| Vacuum Range | From 50 to 550 mmHg |
| Flow (free flow) | 32 l/min |
| Noise level | 70 dB |
| Electrical connection requirements | 100-240 VAC 47-63 Hz 0,75A max; 12VDC 33W max |
| IP Grade | IP12 |
| Battery | 12V DC Ni-MH |
| Operating time | 30 min. ON, 30 min. OFF |
| Battery recharging time | 5 h |
| Container capacity | 800 ml |
| Patient tube lenght | 1,8 m |
| Operating temp. | From 0°C to +45°C |
| Storage temp. | From -40°C to +70°C |
| Umidity | From 0 to 95% |

Med.Device CE marked according to Dir. 93/42/CEE

AMBUJET P 800 D 20G
SC75500C



Registration number

79723

Rev.0 (31/07/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it