

AMBUJET 800 D 10G

Portable suction device
with 10G wall fastener



Space-saving suction device to be used to remove fluids from the airway or respiratory support system. The device creates a negative pressure that draws fluids through disposable tubing that is connected to the collection container. The base performs as a fastener in accordance to EN 1789 standard.

Specific Characteristics

- Low weight and size
- Easy view of the vacuum set thanks to the led indicators
- Classification ISO High Vacuum/ High Flow
- Suction device protected from bumps and falls thanks to the protective shell
- Mechanism of engagement and release from the base extremely rapid
- Direct connection between canister and vacuum pump

Standard Equipment

- 10g fastener with quick unlocking system
- Battery
- Disposable canister
- Disposable patient tube 1,8 m
- Shoulder strap for easy transport

Optional accessories/Consumables

- | | |
|----------|---|
| SC75012E | Lateral bags for Jet |
| SC75014D | Supply cable for cig.lighter 12vcc |
| SC73017E | 220V Power supply |
| RISC005E | Antibacterial filter |
| SC73016E | Sterile suction catheter 1,8m |
| SC73039E | 800 ml disposable canisters (10pcs set) |

Compliant with

- ISO 10079-1
- EN 1789

Technical data

Jet Dimensions	200x225x285 mm
Jet Dimensions w/fastener	200x250x320 mm
Suction device weight	2,4 ± 0,1 kg
Fastener weight	0,8 ± 0,1 kg
Vacuum Range	From 50 to 550 mmHg
Flow (free flow)	32 l/min
Noise level	70 dB
Electrical connection requirements	100-240 VAC 47-63 Hz 0,75A max; 12VDC 33W max
IP Grade	IP12
Battery	12V DC Ni-MH
Operating time	30 min. ON, 30 min. OFF
Battery recharging time	5 h
Container capacity	800 ml
Patient tube lenght	1,8 m
Operating temp.	From 0°C to +45°C
Storage temp.	From -40°C to +70°C
Umidity	From 0 to 95%

Med.Device CE marked according to Dir. 93/42/CEE

AMBUJET 800 D 10G
SC75750C



Registration number

81604

Rev.0 (31/07/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it