

# RES-Q-VAC

## Manual suction device



Compact size manual suction device, for the removal of fluids from the respiratory tract or the support system to breathing and infectious material from wounds. The aspirator creates a vacuum that extracts the fluid through the disposable catheter connected to the collection vessel. Doesn't require power sources and can be operated by hand or foot.

### Specific features

- Lightweight and compact
- Vessel housing with interlocking
- Autoclavable vessel with overflow valve
- Disposable antibacterial filter to limit the risk of contamination
- Aspirator base with non-slip inserts
- Metal foot pedal to provide a greater strength and reliability.

### Accessories/Consumables

RV71102A	RVX 02 – Vessel autoclavabile 0.38lt
RV71104A	RVX 04 - PVC Carrying case
RV71106C	VP FILTER

### Compliant with

EN ISO 10079-2

### Technical data

Dimension	225x165x150 mm
Max vacuum	300 mmHg
Flow	22 l/min
Length suction tube	1300 mm
Vessel capacity	400 ml
Material	PE, PVC, PC, Fe
Weight	1,4 kg

### Class IIa med.dev. according to dir. 93/42/CEE

RES-Q-VAC RED MANUAL  
ASPIRATOR



### RV71120B

CND classification	Z120105
Registration number	331226

Rev.0 (26/07/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
[spencer.it](http://spencer.it)