


## Hivac

### Motor driven vacuum pump



Membrane vacuum pump lubrication free. The device be installed inside of a suitable pneumatic circuit.

#### Specific Characteristics

Low weight and dimensions  
 2 installation holes for M6 screws  
 Stripped connection cables  
 Equipped with input and output connection hoses

#### Technical data

Dimensions	210x170x130 ± 5 mm
Mounting holes interaxis	75 mm
Noise level	97 dBA
Weight	2,65 ± 0,1 kg
Hose	Ø 9 mm
Max flow	53 L/min
Max vacuum	760 mmhg
Power supply	12Vcc
Power	65 W
Absorption	9,5 A

#### Compliant with directive

2004/108/CE  
 1000/14/CE

VACUUM ELECTROPUMP 60 LT.  
 OX02001B



Distributed by Spencer Italia S.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
[www.spencer.it](http://www.spencer.it)