


Haas 15

Motor driven vacuum pump



Suction pump to be installed inside of a suitable pneumatic circuit.

Specific Characteristics

Extremely low weight and dimensions
2 holes for installation
Stripped connection cables

Technical data

Dimensions	115x85x58 mm
Weight	525 g
Max flow	15 L/min
Max vacuum	550 mmhg
Power supply	12Vcc
Absorption	3A
Mounting holes interaxis	47 mm

HAAS 15 - VACCUM PUMP
15LTS. 12V
OX02315B

Distributed by Spencer Italia S.r.l.

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it