

## R-AID Trolley PRO

### Multipurpose backpack with trolley



Emergency trolley extremely versatile and robust, usable also as a backpack. Is characterized by sturdy wheels to ensure high load resistance and by the possibility to be transported in four different ways thanks to carrying handles, shoulder strap, trolley system and shoulder straps.

#### Specific features

Portability without compromise
Made of non-absorbent PVC
4 external pockets with zip closure
1 internal transparent pocket with zip closure
Padded main compartment to provide maximum protection to the content
Strap bands to fix accessory bags
1 main compartment with straps to fix oxygen tank
Aluminum trolley structure
Adjustable and padded shoulder straps
Lumbar belt for a better distribution of the load on the operator
Padded back to increase the comfort
Reflective inserts for greater visibility

#### Included equipment

M-Series bags	
CB00406A	M1 bag – Suitable for adult intubation
CB00402A	M2 bag – Suitable for paed. intubation
CB00403A	M3 bag – Suitable for medicines
CB00405A	M4 bag – Suitable for ventilation
CB00404A	M5 bag – Suitable for infusion

#### Technical data

Dimensions	420x380x620 mm
Weight <sup>(Without pouches)</sup>	3,3 kg
Volume	51 l
Materials	PVC, Al, Nylon
Loading capacity	30 kg

R-AID TROLLEY PRO MADE IN RED PVC WITH POUCHES  
**CB04171B**

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)