

B-BAG

Double bag resuscitator



Double bag self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific features

The double bag transmits the pressure variations to the operator more effectively
Possibility to connect the bag to sources of oxygen supply
Low storage space required

Technical data

| | |
|--------------------|--------------------------|
| Materials | Rubber, Silicone, PC, PE |
| Bag volume | 280 ml |
| Overpressure valve | 40 cmH ₂ O |
| Weight of the kit | 200 g |

Standard equipment

| | |
|---------|---------------------------------|
| EV50007 | EVX 07 – Non reb. Overpr. valve |
|---------|---------------------------------|

Med. Dev. CE marked according to Dir. 93/42/CEE

B-BAG NEONATAL MANUAL
RESUSCITATOR W/OVERPR.VALV
BG07021B



| | |
|---------------------|-------|
| Registration number | 58904 |
|---------------------|-------|

Optionals

| | |
|----------|---------------------------------------|
| OX20131A | OXX30 anti crushing oxygen tube 10pcs |
|----------|---------------------------------------|

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it