

# B-Life

## Reusable resuscitator



Reusable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

### Specific features

- The material of construction and the surface finish increase grip to the operator's hands
- Intended use CHILD and REUSABLE marked on the device to limit risk of incorrect applications
- Connection standard ISO 15/22 mm
- Possibility to connect the bag to sources of oxygen supply
- Low storage space required
- Steam sterilizable

### Standard equipment

EV50007C	EVX 07 - Overpressure and non-rebreathing valve
----------	---

### Optionals

BG03091A	Reservoir
RM20403B	B-LIFE – Silicone face Mask sz.3
OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator

### Technical data

Materials	Silicone, PC
Bag volume	500 ml
Reservoir volume	2500 ml
Stroke volume range	300 ml with 1 hand 350 ml with 2 hands
Overpressure valve	40 cmH <sub>2</sub> O
Maximum oxygen concentration <sup>(with reservoir)</sup>	98%
Max. re expansion rhythm	100 bpm

MD CE marked according to Dir. 93/42/CEE

B-LIFE - SILICONE PEDIATR.  
W/VALVE SOVR.R-VALVE

**BG03023A**

Registration number



21723

Rev.1 (07/06/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)