

B-Life

Reusable resuscitator



Reusable self inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific Characteristics

The material of construction and the surface finish increase grip to the operator's hands
 Intended use ADULT and REUSABLE marked on the device to limit risk of incorrect applications
 Connection standard ISO 15/22 mm
 Possibility to connect the bag to sources of oxygen supply
 Low storage space required
 Autoclavable ^(except reservoir)
 Latex Free

Technical data

Materials	Silicone, medical grade PVC, PC
Bag volume	1600 ml
Reservoir	2500 ml
Stroke volume range	770 ml with 1 hand 900 ml with 2 hands
Overpressure valve	60 cmH ₂ O
Maximum oxygen concentration	95%
Max. re expansion rhythm	45 bpm
Overall weight	485 ±10 g

Standard equipment

BG03091A	Reservoir
EV50011C	EVX 11 Overpressure and non-rebreathing valve
RM20405A	B-LIFE face mask Sz. 5

MD CE marked according to Dir. 93/42/CEE

B-LIFE - SILICONE ADULT
 W/OVERP. RES. R-VALVE AND
 MASK SIZE 5



BG03027A

Registration number

23763

Optionals

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator

Rev.1 (29/07/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it