

PRO B-Life PVC Resuscitator



Disposable self-inflating manual resuscitator for mechanical ventilation of trauma patients or with breathing difficulties.

Specific features

Connection standard ISO 15/22 mm
Possibility to connect the bag to sources of oxygen supply
Low storage space required
Latex Free

Standard equipment

EV50008C	EVX 08 - PED/NEON DISPOSABLE OVERPR.VALVE N/REBR.
RM10891A	Cuffed Mask Sz.1

Optionals

OX20131A	OXX30 anti-crush tube 10pz
BG03091A	Reservoir bag

Technical data

Materials	PVC, PC
Bag volume	500 ml
Overpressure valve	40 cmH ₂ O
Maximum oxygen concentration ^(With reservoir)	98%
Max. re expansion rhythm	100 bpm

MD CE marked according to Dir. 93/42/CEE

PRO B-LIFE - PVC PAED. RES.
W/OVERP. CUFFED MASK

BG02043A



Rev.0 (21/06/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it