

CO B-Life

PVC Resuscitator



Disposable self inflating manual resuscitator kit for mechanical ventilation of trauma patients or with breathing difficulties.

Specific Characteristics

The material of construction and the surface finish increase grip to the operator's hands
 Intended use PEDIATRIC marked on the device to limit risk of incorrect applications
 Connection standard ISO 15/22 mm
 Low storage space required
 Latex Free

Technical data

Materials	Medical grade PVC, PC
Bag Volume	500 ml
Reservoir volume	2500 ml
Overpressure valve	40 cmH ₂ O
Maximum oxygen concentration	95%
Max. re expansion rhythm	100 bpm

Standard equipment

BG03091A	Sacca Reservoir
EV50008C	EVX08 Overpressure/non-rebreathing valve
RM10890A	Cuffed mask Size 1
RM10891A	Cuffed mask Size 2
OX20130A	OXX30 anti crush oxygen tube

MD CE marked according to Dir. 93/42/CEE

CO B-LIFE PEDIATRIC W/VALVE
 CUFF. S1-2 ANTI-CRUSH
BG08005A



Optionals

OX20131A	OXX30 anti crush oxygen tubes 10pcs
BG02001A	Laxan transparent bag for manual resuscitator

Rev.1 (09/01/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it