


MARK PRO

Professional emergency backpack



Professional backpack characterized by extreme versatility thanks to the arrangement of the internal parts, the protection offered to the content and the dual possibility of transport.

Specific features

Equipped with two handles to allow transport also as a bag. The handles are on the upper and lateral sides

1 external pocket with zip closure on the lower part and 1 frontal pocket with access from the side

2 outside pockets with zip closure for telephone or mobile radio

3 transparent pockets with zip closure on the internal side of the main compartment of the backpack

The internal separator has two big pockets with zip closure. One of them is transparent

2 transparent pockets with zip closure on the internal side of the front flap

Aluminum reinforcement on the back

The main compartment is padded to increase protection to the content

The bottom of the backpack is PVC coated to limit the accumulation of dirt and to make cleaning easier

Padded shoulder straps and adjustable in height, with chest and abdominal belts equipped with quick release buckle for a better distribution of loads

Backrest with paddings in the support areas and covered with a mesh to increase transpiration

Reflective inserts to increase visibility

Standard equipmnt

| | |
|----------|--------------------------------|
| CB00406A | M1 SERIES - SOFT PHIAL HOLDERS |
| CB00402A | M2 SERIES - SOFT PHIAL HOLDERS |
| CB00403A | M3 SERIES - SOFT PHIAL HOLDERS |
| CB00404A | M5 SERIES - SOFT PHIAL HOLDERS |
| CB00405A | M4 SERIES - SOFT PHIAL HOLDERS |

Technical data

| | |
|------------------|-------------------|
| Dimensions | 360x520x300±10 mm |
| Backpack weight | 3,13 ± 0,03 kg |
| Overall weight | 4,85 ± 0,03 kg |
| Volume | 38 ± 1 l |
| Loading capacity | 30 kg |
| Materials | Nylon, PE, PVC |

MARK PRO - PROF. BACKPACK RED/BLACK
W/5 POUCHES
CB00465B

Rev.0 (29/06/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it