

## FLEX MASK

Mouth to mouth mask



Disposable face mask with pre inflated pillow, to be used as a protection between the rescuer and the patient during CPR procedures.

### Specific features

Transparent shell to improve patient monitoring  
 Inlet with cap to allow connection of an oxygen supply  
 Hard case for the best protection

### Technical data

Mask dimension	122x95x100 ± 2 mm
Box dimensions	135x110x50 ± 2 mm
Materials	PVC, Nylon
Mask weight	45 ± 5 g
Overall weight	75 g

Ass. Med. Dev. According to dir. 93/42/CEE

FLEX MASK SET WITH  
 ORANGE PVC CASE INS. ON  
 THE BAC  
 TA09020A



Rev.0 (28/03/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)