

FLEX MASK

Mouth to mouth mask



Disposable face mask with pre inflated pillow, to be used as a protection between the rescuer and the patient during CPR procedures.

Specific features

Transparent shell to improve patient monitoring
 Inlet with cap to allow connection of an oxygen supply
 Hard case for the best protection

Technical data

Mask dimension	122x95x100 ± 2 mm
Box dimensions	135x110x50 ± 2 mm
Materials	PVC, Nylon
Mask weight	45 ± 5 g
Overall weight	75 g

Ass. Med. Dev. According to dir. 93/42/CEE

FLEX MASK SET WITH
 ORANGE PVC CASE INS. ON
 THE BACK
TA09019A



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it