

SUMMIT

Multi-purpose emergency backpack/bag



Multi-purpose bag characterized by large volume with small dimensions. It comes with removable pockets anchored by strap.

Specific Characteristics

Main compartment characterized by:

- Detachable semirigid PVC coated divisor with pockets of various sizes and elastic bands to store accessories
- Mesh pockets on the front internal side
- 3 bags of different colors with transparent top to allow an immediate identification of the content
- Adjustable strap keeps the backpack in open position when it is hanged or worn

Secondary compartment with zip closure

Frontal pocket with zip closure

Pocket on the back side with strap closure

Bottom with PVC coating to protect from dirt if the backpack comes in contact with the ground.

Padded shoulder straps with quick-release buckle

Frontal transparent pocket to insert a customized logo

Technical data

Dimensions	470x360x200 ± 30mm
Logo pocket dimensions	265x120 ± 30mm
Weight	1,97 ± 0,1 kg
Volume	25 ± 1 L
Material	Nylon, PVC

SUMMIT ORANGE-MULTI-PURPOSE EMERGENCY
BACKPACK/BAG
CB04028A

Rev.1 (14/11/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it