

## LIENZ KRONE

Professional aluminum case  
for advanced resuscitation equipment



Hard shell case for housing and transport devices generally used in resuscitation procedures.

### Specific Characteristics

|  |
|--|
| Extremely robust metal structure   |
| Wide handle for easy transport   |
| Metal support feet   |
| Rounded corner protectors made of chromed steel  |
| Rapid closure system with lock   |
| 180° opening to maximize access to content   |
| Mobile internal divisors   |
| Thermo-formed compartments to store accessories, one of which is foldable                          |
| Dedicated compartments for oxygen tank and manual resuscitator                                     |
| The transparent separator is fixed by magnet, allows easy view of the content and keeps it ordered |
| 2 strap belt to fix an oxygen tank (max 2 L without ring nut)                                      |
| 3 rings for the attachment of the supplied straps allow to carry the case as a backpack            |
| Equipped with phial holder   |

### Technical data

|                     |                    |
|---------------------|--------------------|
| External dimensions | 530x385x230 ± 5 mm |
| Materials           | AL, Steel, PE      |
| Weight              | 6,5 ± 0,1 kg       |

LIENZ KRONE - ALUMINIUM CASE  
EMPTY  
**CB06002A**

Rev.0 (22/10/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)