

## L-Handle

### Adult traditional light laryngoscope handle



Adult handle to be used with traditional light blades.

#### Specific features

Non-slip surface treatment  
Needs 2 batteries type C not included

#### Blades connection standard

UNI EN ISO 7376

#### Technical data

Materials	Chromed brass
Handle diameter	30 mm
Overall length	150 mm
Handle length	120 mm
Weight (without batteries)	195 g

Class I Med.Dev. according to Dir. 93/42/CEE

L-HANDLE CONVENTIONAL  
LIGHT ADULT  
DG03100A



CND Classification	Z12021003
Registration number	141647

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)