

L-HANDLE

Adult disposable fibre optic laryngoscope handle



Disposable adult laryngoscope handle to be user with fibre optic blades.

Specific features

Non-slip finish

Technical data

Material	PC
Diameter	30 mm

Blades connection standard

UNI EN ISO 7376

DM di Classe I in accordo alla Dir. 93/42/CEE

L-HANDLE DISPOSABLE FIBRE-
OPTIC ADULT

DG03800C

Classificazione CND

Identificativo registrazione



Z12021003

141676

Rev.0 (02/02/2018)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it